FUJIFILM



FUJI COMPUTED RADIOGRAPHY

Extremely clear and optimal imaging founded on extensive experience. Introducing the new epochal and compact FCR from Fujifilm.





FUJIFILM FCR CAPSULA X Specifications

Standard Components:

• FCR CAPSULA X Image Reader (Model: CR-IR 357)

- · AC power cord
- Other System Components:
- · IP cassette type CC, LC
- CR Console
- Image Recorder : DRYPIX 1000/3000/4000/7000

Supplies:

Imaging Plate:

- ST-VI: 14" × 17", 14" × 14", 10" × 12", 8" × 10", 15 × 30 cm (Inch settings) 35 × 43 cm, 35 × 35 cm, 24 × 30 cm, 18 × 24 cm, 15 × 30 cm (Metric settings) IP Cassette:
- Type CC: 14" x 17" (35.4 x 43 cm), 14" x 14" (35.4 x 35.4 cm), 10" x 12", 8" x 10", 24×30 cm, 18×24 cm, 15×30 cm
- Type LC: 35.4 x 124.5 cm, 35.4 x 101.7 cm, 35.4 x 83 cm, 25.2 x 58 cm,
 - 24 x 57 cm

Time Required for IP Feed/Load:

ІР Туре	Required Time
14" × 17" (35 × 43 cm)	Approx. 82 sec.
14" × 14" (35 × 35 cm)	Approx. 72 sec.
10" × 12"	Approx. 66 sec.
8" × 10"	Approx. 53 sec.
24 × 30 cm	Approx. 66 sec.
18 × 24 cm	Approx. 50 sec.
15 × 30 cm	Approx. 65 sec.

The time required for changing the cassette is assumed to be 0 (zero) seconds.
The time required for image erasure is assumed to be for exposure performed at 25 mR of X-ray dose.

Processing Capacity

ІР Туре	When connected to DRYPIX 7000/CR Console Plus
14" × 17" (35 × 43 cm)	Approx. 43 IPs/hr.
14" × 14" (35 × 35 cm)	Approx. 50 IPs/hr.
10" × 12"	Approx. 54 IPs/hr.
8" × 10"	Approx. 67 IPs/hr.
24 × 30 cm (ST)	Approx. 54 IPs/hr.
18 × 24 cm (ST)	Approx. 72 IPs/hr.
15 × 30 cm (ST)	Approx. 55 IPs/hr.

• The time required for changing the cassette is assumed to be 0 (zero) seconds.

The time required for image erasure is assumed to be for exposure performed at 25 mR of X-ray dose.

Time to Print on DRYPIX 7000 through network via CR Console:

- Approx. 135 sec. (35 x 43 cm)
- Time To Display on CR Console:
- 35 x 43 cm: Approx. 52 sec.
- 18 x 24 cm: Approx. 31 sec.

Number of Stackers: 1

Reading Gray Scale: 12 bits

Network: 10 Base T/100 Base TX

Dimensions (W x D x H): 590 x 380 x 810 mm (23" x 15" x 32")


Weight: 99 kg (218 lbs.)

Power Supply Conditions:

Single phase 50-60Hz 120-240V ±10% 5A (max)

Environmental Conditions:

- · Operating Conditions: Temperature: 15-30°C
- Humidity: 15-80%RH (No dew condensation) Non-operating Conditions:
- Temperature: 0-45°C
- Humidity: 10-90%RH (No dew condensation)



IP Cassette Type CC







FUJ!FILM FUJIFILM Corporation 26-30, NISHIAZABU 2-CHOME, MINATO-KU, TOKYO 106-8620, JAPAN

Specifications and PC requirements are subject to change without notice. All brand names or trademarks are the property of their respective owners.